

Activity measure


Measurement of weekly activity

Properties

- Type: Medico
- Groupable: No
- Show in lists: Yes
- Searchable: Yes
- Requires configuration: No


User interface

Field in show mode



Motion 1 gange/uge 2t 15m (135)

Field in edit mode



Motion 1 gange/uge 2 timer 15 minutter

Configuration



Edit field configuration

Disable frequency

Frequency name

CANCEL

[Description of configs]

Developer info

- FeltTypeID: 643
- SQL datatype: DECIMAL(10,4)
- Class name: FieldTrainingSurvey

Revision #2

Created 1 April 2025 14:08:29 by Theis Villumsen

Updated 28 July 2025 14:24:57 by Theis Villumsen